





## Flammability Tester (UL-94)



### Model No. PFT - UL-94

**Flammability Tester** UL94 intends this standard to serve as a preliminary indication of a plastics acceptability for use as part of a device or appliance with respect to its flammability. It is not intended to reflect the hazards of a material under actual fire conditions.

UL 94 flammability testing is the first step toward obtaining a plastic recognition and subsequent listing in the "Plastics Recognized Component Directory" (former known as "Yellow Cards"). UL 94 contains the following tests: 94HB, 94V, 94VTM, 94-5V, 94HBF, 94HF and Radiant Panel.

The following is a brief description of three tests that apply to stock shape products (sheet, rod, tube and film) and how the rating system works. This is not meant as a procedure for running the tests nor as a way of determining the acceptability of a material for a particular application.





# Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India  
 P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com



- Faridabad
- Sonipat
- Kolkata
- Mumbai
- Pune
- Ahmedabad
- Chennai
- Bangalore
- Hyderabad

[www.prestogroup.com](http://www.prestogroup.com)

## Key Specifications:

Display	Digital 7 Segment
Distance between Burner & Specimen	10mm ± 1mm
Specimen Size	Horizontal Burning Test - 1/2"X5" Vertical Burning Test - 1/2"X5" Vertical Burning Test Thin Material 8"X2"
Accuracy	± 2% at full range.
Power	220V, Single phase, 50 Hz
Specimen Holder Angle	45° & 90°
Digital Timer Range	999.9 Second
Digital Timer Least Count	0.1 second
Igniting Flame Length	20mm ± 1

## Highlights:

Flame Auto Stop	Yes
Material	Stainless Steel
Dust Collector Tray	Yes
Cylinder Regulator with Pipe	Yes
Specimen Holder	1 set
LPG Empty Cylinder	Yes
Wire Gauge	Yes
Timer Auto Start	Thought Limit Switch

\* LPG Gas not a part of supply

## Optional

Flammability Tester with Various Standard ASTM-D1230, ISO 3795 & FMVSS 302 also available on request

*Thank you customers for choosing us as your partners in growth !*



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

## OUR OTHER PRODUCTS :-

