





 $\label{lem:condition} Variable \ temperature \ range \ +2... +15 ^{\circ}C \ and \ -15... -25 ^{\circ}C$ Efficient and environmentally-friendly refrigerants

Forced-air or static cooling

Friendly-digital temperature displays with integrated data memory White coated sheet steel housing

LED interior light

Structural specifications

Internal and external structure data	
External material	white
Door material	steel
Internal material	commercial grade moulded
Shelf grids material	wire shelves plastic-coated
Controller data	
Temperature display	external digital
Key lock data	
Key lock	
Alarm security data	
Alarm type	Visual and audible temperature and open door

Model code variations

Interior volume	Cooling system	Temperature refrigerator	Temperature freezer	Description	Model code
350	Static		350 liters, white external material with static cooling	KC350	
	Ventilated		-	350 liters, white external material with ventilated cooling	KC350V
100	Static	+2+15°C	-1525°C =	400 liters, white external material with static cooling	KC400
400	Ventilated		_	400 liters, white external material with ventilated cooling	KC400V



Technical specifications

	KC350	KC400				
Model	KC350V	KC400V				
Temperature data						
Temperature range fridge compartment [°C]	+2+15	+2+15				
Temperature range freezer compartment [°C]	-1525	-1525				
Measures						
Net capacity fridge compartment [liters]	221	269				
Net capacity freezer compartment [liters]	87	87				
External dimensions [WxDxH mm]	600x625x1820	600x625x2010				
Usable shelf area [WxD mm]	440x420	440x420				
Controller data						
Type of control	electronic	electronic				
Cooling data						
Cooling system fridge compartment	static or forced air	static or forced air				
Cooling system freezer compartment	static	static				
Defrost fridge compartment	automatic	automatic				
Defrost freezer compartment	manual	manual				
Illumination data						
Interior light	LED light					
Shelf grid data						
Adjistable shelves	4	5				
Drawers	2	2				
Shelf loading [Kg]	45	45				
Alarm data						
Warning alarm	Visual and audible	Visual and audible				
Sound power data						
Sound power level [dB]	39	39				
Electrical data						
Rated Voltage [V]	220/240	220/240				
Power Frequency [Hz]	50	50				