

The new PREMIUM – ECO HPLC columns from VDS optilab!


VDSpher silica gel with a pore size of 100 Å and a carbon content of 17 % for a very wide range of applications!

This high-quality separation phase is only available in a particle size of 5 µm and in the two most commonly used column dimensions, for which the manufacturing process has been customised and optimised.

Separation columns of the highest quality at amazingly low prices!

P/N VEE181J046250	VDSpher P-ECO 100 C18	5 µm	250 × 4.6 mm	HPLC Column
P/N VEE181J046150	VDSpher P-ECO 100 C18	5 µm	150 × 4.6 mm	HPLC Column
P/N VEE181J046010V	VDSpher P-ECO 100 C18	5 µm	10 × 4.6 mm	Guard Cartridges 5/pk

THE ALL-ROUNDER HPLC COLUMN FOR EVERY LABORATORY!



Individual peak shapes indicate residual silanols, more details in our VDSpher info booklet „Use a better Column“, p. 10